


Link do produktu: <https://www.extramed.pl/czujnik-spo2-typu-nellcor-jednopacjentowy-3m-przylepny-dla-doroslych-wtyk-79-pin-p-927.html>



Czujnik SpO2 typu Nellcor, jednopacjentowy, 3M przylepny, dla dorosłych, wtyk 7/9 pin

Cena brutto	40,50 zł
Cena netto	37,50 zł
Czas wysyłki	24H. W razie braku w magazynie dostawa 2-3 tygodnie (opóźnienia przez COVID-19)
Numer katalogowy	01J-21-01

Opis produktu

Lista referencyjna:	PRODUCENT	Numer OEM
	Covidien > Nellcor	D25
	EnviteC	1002341 (DA-2411-1)
	GE Healthcare > Marquette	2023217-001
	Philips	M1131A , 989803128531
	Unimed	U503-01

Kompatybilność:	PRODUCENT	MODEL
	Artema S&W	Diascope
	Bionet	BM3, BM3 Plus
	Bionics	BPM-200
	Covidien > Nellcor	N-20, N-3000, N-390, N-395, N-6000, N100 (req. preamp cable), N180, N200 (req. preamp cable), NPB-290, NPB-295, NPB-40, NPB-75
	Drager	Infinity Vista, Narkomed, Narkomed 4, Narkomed 6400, SC 6000, SC 7000, SC 8000, SC 9000XL, Vitalert 1000
	Edan	H100B
	GE Healthcare > Corometrics	118, 120, 129, 188, 510, 556
	GE Healthcare > Critikon > Dinamap	1000, 3000, 9710, 9720, Compact, Dash Series, Dinamap Select, Eagle Series, Eagle SMK, Hellige Eagle, MPS, Oxyshuttle, Plus, Pro 300, Pro 400, SMU, SMV, SOLAR Series, TRAM Modular x50SL Series
	GE Healthcare > Marquette	Cardioserv, Dash 2500, Dash 3000, Dash 4000, Dash 5000, Eagle 1000, Eagle 4000,

	i/9500, Procare B40, Solar 8000, Solar 8000M
Goldway	UT4000Apro, UT4000C, UT4000E, UT4000F, UT4000F Pro, UT6000A, Vet 420A
Hunleigh	SC1000
Infinium	Infitron II
Invivo	4500 Plus, 4500 Plus 3, Escort 100 OPT11A, Escort 300A, Escort II OPT30, M12, M12 (Nellcor SpO2), Millennia, Millennia 4500 Plus 2, Omega 1445, Omni-Trak TVS, Prism (Nellcor SpO2), Sat Scout
Ivy Biomedical	405A, 405D, 405P, 405T, Vital-Guard 450 CN
Lumeon	Pulsoksymetry podręczne
Medtronic > Physio Control	Lifepak 12
Mennen	Horizon XL (Nellcor SpO2), MR 1300 Mercury 1200
Mindray > Datascope	Accutorr Plus, Accutorr V, DPM Central Station, DPM1, DPM2, DPM3, DPM4, DPM5, DPM6, DPM7, Expert, Passport, Passport V, PM-50, PM-60, PM-60Vet, VS 800, VS 900
Nihon Kohden	BSM-1102 Life Scope EC, BSM-2304 iPro (Life Scope i), BSM-4102 Procyon LT (Life Scope P), BSM-4104 Procyon (Life Scope P), BSM-4112 Procyon LT (Life Scope P), BSM-4114 Procyon (Life Scope P), BSM-9510 Life Scope M (Nellcor SpO2), OGS-2000 PocketCare (Nellcor SpO2)
Pace Tech	Vitalmax 800 Series (black Nellcor SpO2)
Philips	78352C, 78354C, 78834C, 862108, 862231, 862439, 862474, 862478, 863063, 863069, 863073, 863077, 863266, C1, M1020A, M1020B, M1025B, M1205A, M1350B, M1350C, M1722A, M1722B, M1732A, M1732B, M2475B, M2600A, M2601A, M2703A, M2704A, M2705A, M3000A, M3001A, M3002A, M3500B, M3535A, M3536A, M3926A, M3927A, M3928A, M3929A, M4735A, M8102A, M8105A, M8105AS, M8105AT, MCMS, Merlin Systems, V24, V26, Viridia A1, Viridia A3
Siemens	Sirecust 630, Sirecust 700, Sirem SpO2 Module
Spacelabs	90351-0/6, 90465 (Nellcor SpO2), 90466 (Nellcor SpO2), 90467 (Nellcor SpO2), 90489 (Nellcor SpO2), 90496 / 90369 / 90367 Ultraview (Nellcor SpO2), 90651A-08, IM77

Suntech	wszystkie modele z NIBP 247B włącznie
Welch Allyn	Atlas (Nellcor SpO2), Micropaq (Nellcor SpO2), Propaq 106-EL, Propaq 202-EL, Propaq 204-EL, Propaq 206-EL, Propaq CS (Nellcor SpO2), Spot Vital Signs, Vital Signs (nie-Oxi-Max)

Specyfikacja:

Waga pacjenta	Dorośli > 30 kg
Materiał mocowania	taśma 3M (beżowy plaster)
Materiał kabla	PVC
Wtyk	DB9 7 lub 9pin
Technologia	Nellcor
Certyfikaty	CE, ISO EN 13485:2012, 13485:2012/AC:2012, FDA
Gwarancja	Gdy wada zostanie stwierdzona przed użyciem - wymieniamy na nowe