

Link do produktu: <https://www.extramed.pl/kabel-kompletny-ekg-aami-3-odprowadzenia-klamra-wtyk-6-pin-p-815.html>

Kabel kompletny EKG AAMI, 3 odprowadzenia, klamra, wtyk 6 pin

Cena brutto	216,00 zł
Cena netto	200,00 zł
Czas wysyłki	24H. W razie braku w magazynie dostawa 2-3 tygodnie (opóźnienia przez COVID-19)
Numer katalogowy	2340CI

Opis produktu

Lista referencyjna:	PRODUCENT	Numer OEM
	Criticare	1073/P
	Medtronic > Physio Control	11110-000167
	Unimed	2340P-I

Kompatybilność:	PRODUCENT	MODEL
	AAMI	wszystkie
	ATL	800, 800Plus, Apogee, CX ATL Ultrasound, CX200, CX400, HDI 9/200, Interspec Cardiology XL, Mark IV, Mark VII, Ultramark 4/5/7/9
	Abbott	EX-50, Lifecare 5100, Lifecare 5200
	Advanced Instrumentations	PM-1000, PM-1000EX
	American Optical	26125, 33167, 5260, 5330, Pulsar 4
	Carewell	CPM-8000
	Contec	CMS 9000
	Covidien > Nellcor	C-LOCK, N-1000, N-200, N-250, NPB-3900, NPB-4000
	Criticare	1100, 2200, 504US, 504USP, 506, 507, 507E, 507ELC2, 507S, 508, 602-4B, 602-4USP, 8100E, 8100EP, nCompass, nGenuity 8100E, Poet 2, Scholar
	Datex Ohmeda	enGuard CM5
	Digicare	LW600
	Drager > Air Shields	HRRM-71-1, HRRM-71-2, Porta-Fib I, Porta-Fib II, Porta-Fib III, System 5 (Nowy), System 6 (Noworodki)
	Drager	Vitalert 2000, Vitalert 3000
	Edan	M50, M8, M80, M8A, M8B, M9, M9A
	GE Healthcare > Corometrics	556

GE Healthcare > Critikon > Dinamap	8100T, 8700/9700, MPS Select, Pro 1000
Goldway	G40, UT4000A, UT4000Apro, UT4000B, UT4000C, UT4000F, UT4000F Pro, UT4800, UT6000A, UT6000F
Guthan	1000
Infinium	Infitron II, Omni, Omni II
Invivo	101, 102, 200, 300, 360SL, 450SL, 521A, Escort 100, PortaPak 80, PortaPak 90/1000, Prism
Ivy Biomedical	101 (sprzed 1989), 402, 700
LSI	261, LS5
Mediana	YM6000
Medtronic > Physio Control	1440/R, Lifepak 10, Lifepak 5, Lifepak 6, Lifepak 6S, Lifepak 7, Lifepak 8, Lifepak 9, VSM-1ESF, VSM-2, VSM-2ESF, VSM-3
Midmark > Cardell	9403V, 9405, Max-1, MAX-12 Duo HD
Mindray > Datascope	MEC 1000, MEC 1200, MEC 200, MEC 2000, PM 7000, PM 8000, PM 9000
Nihon Kohden	7102A, BSM-7102A, BSM-7300A, BSM-8108A, BSM-8301A, BSM-8302A, BSM-8500A, BSM-8502A, BSM-8800A, Life Scope 12, Life Scope 6, Life Scope 7, Life Scope 8, Life Scope 9, MU-831RA, OEC-6102A (s/n>20156), OEC-6105A (s/n>20226), OEC-7102A, OEC-8108A, TEC 5531E, TEC-5200A (s/n>20421), TEC-7100A, TEC-7200A
Novamatrix	902, 903, 910, Life Defense 1
Omron > Colin	BP308, BP408, CBM3000, M20, M30, MS20
PPG	Guardian 1001, N1001
PaceTech	2100, 2200, 4000, 4100, MiniMax 4000CL, Minipack 300, Minipack 3000/3100, Modular/4100CL, Vitalmax 4000CL, Vitalmax 800 Plus
Philips	iE33 xMatrix Ultrasound CX50, M3921A A1, M3923A A1, M3925A A1, M3927A A3, M3929A A3
Smiths Medical > BCI	3101, 3404, 6100, 6200, 9100, 9200, Advisor, Autocorr: Surgivet
Spacelabs	67112, 90300, 90303, 90303B, 90400, 90600, 90603A, 90700, Medic 3, Medic 4, Medic 5, Medic 6, PCMS, Series 500, Series 600, Series 900, TEK-208/408/511/512/514/521
Tektronix	4045 Defib/208/400/500 Serie
Welch Allyn	Propaq 206 EMPSC
Zoll	D900, E Serie, M Serie,

Specyfikacja:		PD1200, PD1400, PD1600, PD1700, PD2000, R Serie
	Złącze proksymalne	Klamra
	Długość całk. kabla	2.80 m
	Materiał kabla	TPU
	Wtyk	6 pin
	Kod kolorów	IEC
	Certyfikaty	CE, ISO EN 13485:2012, 13485:2012/AC:2012, FDA
	Gwarancja	6 miesięcy