

Link do produktu: <https://www.extramed.pl/lista-referencyjna-producentnumer-oem-unimedu203s-01-p-39.html>

Czujnik SpO2 typu Nellcor, silikonowy dziecięcy, wtyk 7 pin, kabel 1 m

Cena brutto	324,00 zł
Cena netto	300,00 zł
Czas wysyłki	24H. W razie braku w magazynie dostawa 2-3 tygodnie (opóźnienia przez COVID-19)
Numer katalogowy	01K-310S-01
Kod producenta	CR001-5301HS

Opis produktu

Lista referencyjna:	PRODUCENT	Numer OEM
	Unimed	U203S-01

Kompatybilność:	PRODUCENT	MODEL
	Artema S&W	Diascope
	Bionet	BM3
	Bionics	BPM-200
	Covidien > Nellcor	N-20, N-3000, N-390, N-395, N-6000, N100 (wymagają kabla przedwzmacniacza), N180, N200 (wymagają kabla przedwzmacniacza), NPB-290, NPB-295, NPB-40, NPB-75
	Drager	Infinity Vista, Narkomed, Narkomed 6400, SC 6000, SC 7000, SC 8000, SC 9000XL, Vitalert 1000
	Edan	H100B
	GE Healthcare > Corometrix	118, 120, 129, 188, 510, 556
	GE Healthcare > Critikon > Dinamap	1000, 3000, 9710, 9720, Compact, Dash Series, Dinamap Select, Eagle Series, Eagle SMK, Hellige Eagle, MPS, Oxyshuttle, Plus, Pro 300/400, SMU, SMV, SOLAR Serie, TRAM Modular x50SL Serie
	GE Healthcare > Marquette	9710, 9720, Compact, Dash Series, Dinamap Select, Eagle Series, Eagle SMK, MPS, Oxyshuttle, Plus, Pro 300/400 (sin-Nellcor OxiMax), SMU, SMV, SOLAR Series, TRAM Modular x50SL Series
	GE Healthcare > Marquette	Dash 2500, Dash 3000, Dash

	4000, Dash 5000, i/9500, Procure B40, Solar 8000, Solar 8000M
Huntleigh	SC1000
Infinium	Infitron II
Invivo	Escort 100 OPT11A, Escort 300A, Escort II OPT30, M12 (Nellcor SpO2), Millennia 4500 Plus 2, Prism (Nellcor SpO2)
Ivy Biomedical	405A, 405D, 405P, 405T, Vital-Guard 450 CN
Lumeon	Pulsoksymetry podręczne
Medtronic > Physio Control	Lifepak 12
Mennen	Envoy, Horizon XL (Nellcor SpO2), Mercury 1200, MR 1300
Mindray > Datascope	Accutorr Plus, Accutorr V, DPM Central Station, DPM1, DPM2, DPM3, DPM4, DPM5, DPM6, DPM7, Expert, Passport, Passport V, PM-50, PM-60, PM-60Vet, VS 800, VS 900
Nihon Kohden	BSM-1102 Life Scope EC, BSM-2304 iPro (Life Scope i), BSM-4102 Procyon LT (Life Scope P), BSM-4104 Procyon (Life Scope P), BSM-4112 Procyon LT (Life Scope P), BSM-4114 Procyon (Life Scope P), BSM-9510 Life Scope M (Nellcor SpO2), OGS-2000 PocketCare (Nellcor SpO2)
PaceTech	Vitalmax 800 Series (czarny Nellcor SpO2)
Philips	78352C, 78354C, 78834C, 862108, 862231, 862439, 862474, 862478, 863063, 863069, 863073, 863077, 863266, C1, M1020A, M1020B, M1025B, M1205A, M1350B, M1350C, M1722A/B, M1732A/B, M2475B, M2600A, M2601A, M2703A, M2704A, M2705A, M3000A, M3001A, M3002A, M3500B, M3535A, M3536A, M3926A, M3927A, M3928A, M3929A, M4735A, M8102A, M8105A, M8105AS, M8105AT, MCMS, Merlin Systems, V24, V26, Viridia A1, Viridia A3
Siemens	Sirecust 630, Sirecust 700, Sirem SpO2 Modul
Spacelabs	90351-0/6, 90465 (Nellcor SpO2), 90466 (Nellcor SpO2), 90467 (Nellcor SpO2), 90489 (Nellcor SpO2), 90496 / 90369 / 90367 Ultraview (Nellcor SpO2), 90651A-08, IM77
Suntech	wszystkie modele włącznie z NIBP 247B
Welch Allyn	Atlas (Nellcor SpO2), Micropaq (Nellcor SpO2),

Propaq 106-EL, Propaq 202-EL, Propaq 204-EL, Propaq 206-EL, Propaq CS (Nellcor SpO2), Spot Vital Signs, Vital Signs (nie-OxiMax)
--

Specyfikacja:

Waga pacjenta	Dziecięcy 3-15 kg
Długość kabla	1.0 m
Materiał kabla	TPU
Wtyk	DB9 7pin
Technologia	Nellcor OxiSmart > Covidien
Certyfikaty	CE, ISO EN 13485:2012, 13485:2012/AC:2012, FDA
Gwarancja	12 miesięcy